


1. Measurement & Error

purpose. measure Density.

$$\text{Density} = \frac{\text{mass}}{\text{volume}}$$

unit. g/cm^3


★ Volume of Shapes:



$$V = \pi r^2 h$$



$$\frac{4}{3} \pi r^3$$



$$L \times W \times \text{thickness}$$


Apparatus

1. Triple beam Balance



to find mass

2. Vernier Calliper



to find r or L or W

$$\frac{\text{Zero value} + \text{line value} \times 0.01}{10} \text{ from bottom}$$



$$\rightarrow \frac{25 + 5 \times 0.01}{10} \rightarrow \star \text{ the scale is 10 points it's mm}$$

→ to cm

or just divide all by 1000
mm → m

3. Micrometer Calliper



to find thickness



$$\rightarrow \frac{6.5 + 8 \times 0.01}{10}$$

$$\frac{\text{the not covered value} + \text{Zero value} \times 0.01}{10}$$

10

procedure

1. find Mass using T.B.B.
2. Measure Volume using Vernier Caliper / micrometer
3. Calculate Density using calc.
4. Calculate error using calc.

Relative Error

$$\frac{\text{real val} - \text{measured}}{\text{real}} \times 100$$

Conclusion.

measured Density

2. Thermal

purpose: Find Thermal Expansion Coefficient.

$$\alpha = \frac{\Delta L}{\Delta T \cdot L}$$

unit: $\frac{1}{^{\circ}\text{C}}$

Apparatus

1. Metal tube
2. Steam generator and water
3. Ruler
4. Ohm's meter and wires
5. Thermostat / dial gauge

procedure

1. open the Steam generator after adding water
2. measure the length L of the Tube (from middle of Black circle Left to right)
3. measure R from ohm meter then look in the table to find Temperature
4. Apply the Steam then measure resistance from ohm's meter and find final Temperature from table.
5. measure ΔL directly from digital gauge
6. calculate thermal expansion $\frac{\Delta L}{\Delta T \cdot L}$



	DATA				CALCULATIONS		
	L(mm)	R _m (Ω)	ΔL(mm)	R _{hot} (Ω)	T _m (°C)	T _{hot} (°C)	ΔT(°C)
Copper							
Brass	750	10050	1.05	770	25°	96°	71
Aluminum							

ruler → L
 from reading → R_m = reading × 10³ Ω
 from circle directly after steam → R_{hot}
 after the steam → from table
 after the steam from table → change in T
 then search table for T_m = °C
 when water comes from tube take values (after they are stable)

Red → Temp (1st connection)
 Blue → Com (Sec connection)

3. Energy Conservation

purpose. Approve that Energy is conserved

$$KE_i + PE_i = KE_f + PE_f$$

$$KE = \frac{1}{2}mv^2$$

$$PE = mhg$$


Apparatus

1. pasco software
2. motion sensor
3. interface
4. free fall system



procedure

1. pasco settings ² → chose graphs + table
Data Summary
- motion
- velocity
position → settings → N. Format → 3 → Okay
* for position and velocity



2. connect motion sensor with interface

3. Drop the ball on the sensor

4. collect data

- The Table will fill itself
- look for value when $V = 0.00$
- take the value of position
- then calculate $\frac{3}{4}$ and look for the closest value of p in table
- take the value of V
- Same thing to find V when $p = \frac{1}{2}$
- and $p = \frac{1}{4}$



+ same connection
in details in
Free fall exp.

calculate KE and PE.

Conc.

- total Energy at different height are equal
- total Energy in closed system is conserved.


4. Free Fall

purpose. Study the relationship between Time and displacement on freely falling object, measure g .

$$y_i = \frac{1}{2} g t^2$$

Apparatus.

1. Drop box
2. Control box
3. Time switch
4. adapter
5. Time of flight
6. Ball
7. Ruler
8. cable
9. timer




procedure.

1. connection
2. use ruler to measure initial distance y_i from the bottom of the ball to the Time of flight.
3. drop the Ball and read the timer

4. calculate the acceleration of gravity
plug the time when height = 80, 60, 40, 20, 0
and find t^2 .

plug the point on graph



find slope between any two points (not the ones plugged).

The real value of $g = 9.81$
find relative error.

$$g = \frac{2}{\text{slope}}$$

Conc.

calculated acc of gravity.

5. vector

purpose. Finding resultant force F_R .
resultant angle θ_R .
Equilibrant angle θ_e .

$$F = ma$$

$$F_x = F_1 \cos \theta + F_2 \cos \theta$$

$$F_y = F_1 \sin \theta + F_2 \sin \theta$$

$$F_R = \sqrt{(F_x)^2 + (F_y)^2}$$

$$\theta_R = \tan^{-1} \left(\frac{F_y}{F_x} \right)$$

$$\theta_R = \theta_e - 180$$

$$\theta_e = \theta_R + 180$$

Apparatus.

1. Force table
2. mass & hanger. Set
3. protractor
4. Graph paper
5. pencil
6. ruler

Experimental [always coming]

☆ $m_1, \theta_1, m_2, \theta_2$ will be given.

1. add them to the table.
(the blue hanger is 5g)
2. Try different m and θ until the circle is centered
3. now u have m_e, θ_e
4. find F_R

$$|F_R| = |F_e| = m_e g$$

$$\theta_R = \theta_e - 180$$

Graphical

1. use protractor  to draw θ_1 and use ruler to draw m_1 . (the length is the grams $\div 10$)
if $m = 45g$ the length $\div 10 = 4.5 \text{ cm}$. 
start from zero point
2. Same for m_2, θ_2 but start from wherever you end m_1, θ_1
3. connect them to find m_R, θ_R . the θ take it from zero point 
 $m_R = \text{cm} \times 10 \rightarrow g$
 $\theta_e = \theta_R + 180$

Theoretical

apply Formulas

$$F_x = F_1 \cos \theta + F_2 \cos \theta$$

$$F_y = F_1 \sin \theta + F_2 \sin \theta$$

$$F_R = F_E =$$

$$\sqrt{F_x^2 + F_y^2}$$

$$\theta_R = \tan^{-1} \left(\frac{F_y}{F_x} \right)$$

$$\theta_e = \theta_R + 180$$

$$\text{Relative Error:} = \frac{\text{real value (from exp)} - \text{this}}{\text{real}} \times 100$$

Conc.

calculated resultant force for two other forces

6. Projectiles

purpose. determination of the horizontal range (R) of projectile shot for various velocities and angles.

Zero launch angle:

$$t = \sqrt{2 \frac{h}{g}}$$

$$R = v_0 \sqrt{2 \frac{h}{g}}$$

unit: m

General launch angle

$$R = \left(\frac{v_i^2}{g} \right) \sin 2\theta$$


Apparatus.

1. projectile launcher
2. Smart timer
3. 2 photo gate heads
4. metric measurement tape
5. Carbon - white paper
6. Adapter
7. Loading Rod
8. Ball
9. tape + scissors



procedure.

1. connect smart timer with two photogates
2. connect smart timer with adapter
3. put ball inside projectile (1-3 clicks)
4. using ruler measure initial height and d to measure initial velocities and time



5. check angles

draw
line on
Table

6. Turn on Smart timer
7. pull the ball and mark where will it first hit the carbon paper
(first try to see where u should place the white paper.)
8. measure horizontal Range using metric measurement tape.



Relative Error:

$$\frac{\text{real value. (R calculated)} - R_{\text{measured (Experimental)}}}{\text{real}} \times 100$$

- Gate 1. = Short
Gate 2. = Medium
Gate 3. = Long.

For General Launch Angle:

plot a graph between R measured and it's angle
R calculated and it's angle



Conc:

used the equation of motion to determine R of projectile.

7. Newton Law

purpose. Finding mass of a car using newton's 2nd law.
By studying the relationship between a and F of the car

$$F = m \cdot a$$

by force sensor $\swarrow$ $\downarrow$ our goal $a = \frac{v}{t}$

Apparatus.



1. force sensor
2. dynamics system the road of the car
3. computer interface
4. pas car
5. pasco software
6. smart pulley
7. mass and hanger set
8. motion sensor



procedure.

☆ sitting $\rightarrow$ choose 3
+ change sign

1. open pasco - chose 2 graphs + table.




2. connect force/motion sensors to interface.

3. add hanger and masses to the car.



4. read the force and acceleration (slope of v, t)




$\rightarrow$ force (mean)

$$\text{slope } \frac{\Delta F}{\Delta a}$$

$\hookrightarrow$ the experimental m

compare it with
real m. to find
Relative Error.



$\rightarrow$ velocity (same time,
linear slope)

8. Friction Force

purpose. find coefficient of static friction
find coefficient of kinetic friction

$$\mu_k = \frac{f_k}{N}$$

no unit

$$\mu_s = \frac{f_s}{N} \rightarrow \text{Static friction } f$$

$N = W = mg$

Apparatus.

1. force sensor
2. Triple beam balance
3. Pasco software
4. computer interface
5. Friction Accessory + String
6. Bar masses
7. Motion Sensor

procedure.

1. open pasco + choose 2 graphs + Table
2. Connect motion / force sensor to interface
3. connect force sensor with friction Accessory
face the motion sensor to " "
4. read the f_s (Max) from the graph
read the f_k (mean)



Relative error $\rightarrow$ depending on the material see the real value

Conc. $\mu_s > \mu_k$

4. Spring Force

purpose: find Spring Constant k

$$F = k \Delta x \rightarrow \text{displacement}$$

↓
our goal

$$\text{Slope } \frac{\Delta F}{\Delta x} = k$$

Apparatus:

1. Force Sensor
2. Pasco
3. com interface
4. clamp
5. ruler
6. Springs

procedure:


1. open pasco, chose graph + Table Force + d.
 - choose KEEP mode
 - change force sign
 - press zero on sensor

2. connect force sensor with inter face

3. connect the spring with the clamp + the force sensor, fix the ruler

4. Expand the spring 10cm each time (press keep)

find k from the slope.



Long Shiny Doll } Expand

Heave Light } Compress

conc. found k for different kinds of springs

10. momentum

purpose. investigate the conservation
of momentum & kinetic energy
of a system in collision.

Elastic

$KE_i = KE_f$
Kinetic energy is
conserved

Inelastic

K Energy is not
conserved

P is conserved in both

$$P = m \times v$$

$$KE = \frac{1}{2} m v^2$$

$$\star P_{total,i} = P_{total,f}$$

Apparatus.

1. Mass balance
2. Pasco
3. com interface
4. Dynamic track
5. two Pas cars
6. Smart pulley
7. 250 g mass
8. Rotary motion sensors




procedure.

1. open pasco, choose 2 graphs
2. Connect rotary-motion sensors with interface
3. Do the collisions when $m_1 = m_2$
 $m_1 \neq m_2$

Elastic Inelastic

4. From graph take values of $V_i - V_f$

$\star$ you have to add sensor
manually (digital-RM)
1 and 2



moving car

stable car

Highlight + find [min
max

$$P = m \times v$$

add P_i and P_f to find P_{total}

Relative error the initial P is
the real value

Conc. the momentum
in elastic / inelastic is conserved.

☆ if
u got
negativ
values
exchange
yellow
and black
gates.